Allard ToeOFF® Family A Product for Every Need!

Condition		Ypsilon™	ToeOFF®	BlueROCKER™	Custom ToeOFF®	СОМВО™
Footdrop	Mild					
	Complete					
M-L Stability	Stable					
	Unstable					
Body Type	Thin-Average					
	Average-Stocky					
	Stocky-Obese					
Range of Motion	Normal					
	Hyper/Hypo					
Proximal Deficits	None					
	Mild					
	Moderate					
	Severe					
Spasticity	None					
	Mild					
	Moderate					
	Severe					
Muscle Deficit	Balanced					
	Unbalanced					
Partial Foot Amputations	Тое					
	Transmet					
	Lisfranc					
	Chopart					
The DRIGINAL		P			1	
bon composite,						
י Ioor react						

carbon composite, floor reaction, dynamic response **AFO!**

The following are only "beginning" guidelines for product selection based on various patient conditions. Take into consideration different patient aspects like spasticity, stability, balance, need for support, activity, lifestyle, body type and other factors. This should always be done on an individual basis with the fitter taking into consideration the unique situation and needs for the individual patient.

Call US for your free slide chart!

www.allardusa.com

ALLARD USA, INC. 300 Forge Way, Suite 3 Rockaway, NJ 07866-2056 info@allardusa.com Toll Free 888-678-6548 Fax 800-289-0809