



Control hip alignment
Increase levels of function

S.W.A.S.H.[®]

Sitting Walking And Standing Hip orthosis



CALL TODAY TO RECEIVE OUR NEW CLINICAL & FITTING MANUALS

SWASH is only sold to orthotists who have been trained to fit the orthosis.



Original



Low Profile II

Two options for upright selection



More abduction in sitting

Increased abduction offers maximum opportunity for muscle lengthening in the sitting position, and is also recommended when there is need to influence the trunk toward more spinal extension.



Less abduction in sitting

Reduced abduction better accommodates wheelchairs with narrow seat structures and is better tolerated by children with tight adductor muscles.



Your O&P Partner