ORION More leaturing

The ORION rapidly selects support as needed to safely negotiate environmental obstacles encountered during everyday walking. The knee is versatile, adapting for slow to fast walking speeds without being restrictive. This reduces the need for increased muscle input and leads to a reduction in energy expenditure throughout the day.

The OnBoard Smart Programming is intuitive, resulting in finely tuned gait parameters for enhanced symmetry. The option to use the ORION Software Suite also allows the practitioner to adjust and monitor kinetic and kinematic gait parameters with the ability to save and record these settings to a database.

All of these features give you a microprocessor knee that is solid, stable and dependable. The only thing left to do is to get buck house

To learn more, stop by our booth at the 2013 Amputee Coalition National Conference

Visit www.endolite.com for more information



www.spsco.com 800.767.7776

www.endolite.com 800.548.3534 Scan for video

endolite

