

First-Class Liner

The special silicon has a caring effect on the skin due to the substances contained. This guarantees the patient an optimal fit, an ideal level of comfort and a secure base for the prosthesis. The Jersey material is ladder-proof thanks to the weave and can be cut with our scissors (having corrugated cutting edges) at any position you wish.

Area of application: Specially suitable for geriatric patients and for prosthesis users who are either not very active or are moderately active due to ease with which the liner can be rolled up and down

Art No. 1S5 & 1S6
 Sizes 16-18-20-22-23,5-25-26,5-28-30-32-34-36-40-45

Silicon
 100% Compound crosslinked silicon, cares for the skin without using solvents, long-lasting and durable

Contact surface
 optimal grip, high degree of comfort for the wearer


Easy roll cover
 the silicon is given optimal protection through the high-quality jersey cover and ensures that the patient can roll the liner up and down with ease

Direct power transmission
 through the connection of the connection disc and internal laminate

Distal connection
 with metal thread
 * Also available without distal connection

1S5

1S6



Circle # 46 on Readers' Service Card



SPS
 ORTHOTIC PROSTHETIC SUPPLIES
 Ph 800.767.7776
 Fax 800.869.7776

KNIT-RITE®
 120 Osage Avenue
 PO Box 3900
 Kansas City, KS 66103
 Web: www.knitrite.com
 PH: 800.821.3094
 FX: 800.462.4707

PEL
PEL SUPPLY CO.
 4666 Manufacturing Rd.
 Cleveland, OH 44135
 Ph: 800-321-1264
 Fx: 800-222-6176
www.pelsupply.com

For free information contact

Euro International, Inc.
 1.800.378.2480

www.eurointl.com • sales@eurointl.com
 Ph: 1.813.246.5995 • Fax: 1.813.246.5998
 5906 Breckenridge Pkwy. Suite G • Tampa, FL 33610

