

Stance - Control Knee™

Designed for people who have knee instability due to quadriceps weakness or absence. Results in a more normal and energy-efficient gait than traditional locked knee joints.



Post Polio

Upcoming Courses

July 19

Little Rock, AR
Horton's Orthotic Lab

September 18

Chicago, IL
Academy Advanced Training Seminar

October 11

White Plains, NY
New York State Chapter AAOP

October 21

Reno, NV
AOPA National Assembly

SCOKJ™ Facts



Prevents knee collapse regardless of flexion

Stance control can be triggered by weightbearing or joint motion, according to patient needs

Provides three operating modes:
Automatic Stance Control, Free Motion, and Manual Lock

Suitable for metal, thermoplastic, laminated, and hybrid orthoses

Fabrication and adjustments are straightforward and require no expensive special tooling

Can be used with bilateral knee paralysis or paresis

Horton Technology Inc.

Innovations In Orthotics and Prosthetics

5220 West 12th * Little Rock, AR * 72204

1-866-663-3970

www.stancecontrol.com