


SOLUTIONS THAT FIT™

PERFORMANCE SOCKS & SHEATHS

A/K BRIM SHEATH

For above knee patients, for all socket styles. Designed to separate and protect the soft tissue in the top brim of the socket. Used to bridge the gap between the top edge of any roll on liner and the groin. Will work on pin suspension and seal in liners.

Brim sheath uses stretchy, slick nylon material, constructed in two angled sections to allow the fabric to fit the high lateral and lower medial shape appropriately. Includes gripper elastic to secure to outer socket. Patent pending.



A/K Brim Sheath

KNIT-RITE LINER-LINER® PROSTHETIC SOCK WITH X-STATIC®

A prosthetic sock designed to be worn under a suspension liner next to the skin for skin protection.

- Relieves skin shear irritations and improves comfort with liners
- Tested and proven beneficial for amputees who wear liners as the primary interface
- Ultra-thin material minimizes compromise to liner's suction grip on skin
- Moisture-wicking fibers for a cool and comfortable garment
- X-STATIC silver fibers assist in preventing odor in the garment
- Washable interface keeps liner cleaner and helps control skin irritations
- High stretch for excellent fitting qualities



Knit-Rite Liner-Liner®
Prosthetic Sock with
X-STATIC®

Circle # 11 on Readers' Service Card

For more information contact Knit-Rite at 800-821-3094 or e-mail customerservice@knitrite.com.

Also distributed by:

KNIT-RITE[®]
INCORPORATED

CONNECT WITH US ON THE WEB!