



# M-**BRACE**™

soft wear for the body



**KNEE #46**

## KNEED SUPPORT?

### OA KNEE BRACE WITH RANGE OF MOTION

In your quest to live an active lifestyle, the #46 knee brace brings your comfort to a new level with its elegant and light design. Available for the right or left knee, this brace alleviates knee pain by reducing the weight on your medial or lateral compartment by using the three point pressure system. It is composed of two aluminum reinforced semi-rigid straps that surround the distal region (above the knee) and the proximal region (below the knee). These straps are to be custom moulded to allow for an optimal torsion support. Additionally, these straps are bound together by a patented articulated aluminum support that is easily adjustable for flexion and extension and accessible within a zipped pocket.

Approved for L-1843

Use this knee brace for the following indications:

- Mild to moderate medial compartment OA without instability.
- Meniscus procedures.
- Chondral defect procedures.
- Moderate lateral compartment OA without instability; it is possible to reduce the load on the lateral compartment by using a "right" brace on a left leg or a "left" brace on a right leg.
- Available in left or right.

## BRACE YOURSELF

### AIR BACK BRACE

This Lumbar-Sacral Brace acts as an immobilizing device for the lumbar-sacral region and its' primary function is to reduce pain caused by lumbago, sciatica and ischialgia.

Made with semi-rigid fabric, it enjoys a structure made with eight selfmodeling and flexible steel stays on which a number of rigid straps slide. When the locking system is pulled and secured at the front, the inside straps stiffen the entire structure, thus effectively immobilizing the lumbar-sacral region.

Two rigid thermo-formable panels for lumbar and abdominal application are included. These panels produce intra-abdominal pressure to reduce load on discs thus reducing pain and improving post-operative outcome when used with the corset in the appropriate position. The supported region extends from the sacrococcygeal joint to L-1 vertebra.

Approved for L-0637, L-0631, L-0627



Use this back brace for the following indications:

- Lumbago, sciatica and ischialgia.
- Lumbar-Sacral pain of other origin and nature.
- Disc complaints.
- Spondylo-arthritis.
- Osteoporosis related lumbar vertebrae complaints.
- Para-vertebral lumbar muscle spasm.
- Light trauma of the lumbar-sacral column.

### MIH International LLC

Toll Free Tel: 1 866-374-4050  
Toll Free Fax : 1 866-374-4051  
sales@m-brace.com  
www.m-brace.com

MADE IN ITALY



Contact us today for more information.