



## PIPDriver...

Custom device driven by an intact PIP joint with enough residuum to engage the ring.



Custom device driven by an intact MCP joint with enough residuum to engage the ring.





Custom device driven by intact MCP and CMC joints in the thumb. However, patients with more proximal amputations may still benefit from the device as long as their residuum can engage the ring.

## **GripLock** Finger...

Elegant solution for amputations proximal to the MCP joint. Positioned using a surface or their opposing hand, and then return to full extension by releasing the dorsal pawl or moving to the full flexion release point.











