



NAKED
PROSTHETICS™

GripLock Finger™

Naked Prosthetics' functional solution for amputation proximal to the MCP.



Positioned using a surface or their opposing hand, and then return to full extension by releasing the dorsal pawl or moving to the full flexion release point.

npdevices.com

The **GripLock Finger** device can be worn with or alongside our **MCPDriver**, **PIPDriver** and/or **ThumbDriver**, and provides the same colors and design choices our users love. The **GripLock Finger** is used in conjunction with a custom-fabricated socket and to be fit on a patient by a certified prosthetic clinician.



PIPDriver™

Custom device driven by an intact PIP joint with enough residuum to engage the ring.



MCPDriver™

Custom device driven by an intact MCP joint with enough residuum to engage the ring.



ThumbDriver™

Custom device driven by intact MCP and CMC joints in the thumb. However, patients with more proximal amputations may still benefit from the device as long as their residuum can engage the ring.



For more information: (888) 977-6693 | (360) 915-9724 | fax (360) 918-8619 | info@npdevices.com or visit npdevices.com

© 2021 Naked Prosthetics. All Rights Reserved.

