

# Versatility and support...



## ...the perfect combination for lower mobility patients

Put the 1M10 Adjust foot together with the 3R93 knee for an outstanding blend of security and flexibility.

Designed specifically for lower mobility users, the Adjust foot provides a stable center of gravity – plus easy heel stiffness adjustment to match your patient's needs. The versatile 3R93 knee readily adapts as your patient's mobility increases, converting from manual lock mode to a friction brake knee joint with stance flexion.

Your patient gains security, comfort and confidence from this winning combination. Contact your Otto Bock sales representative at 800.328.4058.

[www.ottobockus.com](http://www.ottobockus.com)



**Otto Bock**<sup>®</sup>

QUALITY FOR LIFE