

Manufacturers of The ONLY Multi Podus® System



Hip Abduction Orthosis

Provides stabilization for patients prone to chronic hip dislocation. Flexion/extension ROM control in 10° increments. Abduction-adduction adjustments in 5° increments. Universal right or left.



Resting Hand Orthosis

(left) The Resting Hand Orthosis offers a functional resting hand position following stroke, trauma, injury, as secondary casting post-surgically, and can be used in the treatment of rheumatoid arthritis, arthroplasties, burns, and early wrist and finger contractures.

Dorsal Carpal Tunnel Splint

The Dorsal Carpal Tunnel Splint limits wrist flexion and extension during repetitive hand motion.



Geriatric Hand

The Geriatric Hand is indicated following stroke, arthritis, trauma, or hand injury, for the treatment of finger, wrist, and thumb contractures.

MPO 2000® Active

Treats ankle contractures, eliminates pressure on the heel and toes. Washable, moisture "wicking" foam liner for comfort. Removable, limited ambulation platform. No tools needed.



Post Operative Pin (POP) Knee Brace

The Post Operative Pin Knee Brace is recommended for post operative immobilization and protected range of motion associated with ACL, PCL, MCL, and LCL surgeries. Available in drop lock hinge and heavy duty versions.



Post Operative Pin (POP) Arm Brace

The Post Operative Pin (POP) Arm Brace is designed for post surgery and/or related elbow injuries where limited ROM or locked control of the elbow is desired.



To order, phone
(727) 573-1595 or (800) 627-1595

View the entire RCAI product line at
www.rcai.com

Made in the USA