



## WalkAide® System Myo-Orthotics Technology™

Represents the merging of orthotic technology, which braces a limb, with electrical stimulation, which restores specific muscle function. Moving beyond bracing, Myo-Orthotics Technology devices restore the functionality of an impaired limb by recreating a natural nerve-to-muscle response using electrical stimulation. WalkAide utilizes advanced sensor technology to analyze leg & foot movement. WalkAide then utilizes Myo-Orthotic Technology to send electrical pulses to the peroneal nerve which prompts the muscles in the leg to raise the foot at the appropriate time.

## 2007 WalkAide® Training Schedule

DATE

## Sep. 10-11 Detroit, MI Oct. 15-16 Newport Beach, CA Oct. 22-23 Jacksonville, FL Nov. 5-6 Albany, NY

LOCATION

Come view the New Patient Testimonial Video in Booth #300 at the AOPA National Assembly in Las Vegas.



Experience Our Commitment

Tel: 800.767.7776 x 3 www.spsco.com

Fax: 800.869.7776 estore.spsco.com