



WalkAide® System Myo-Orthotics Technology™

Represents the merging of orthotic technology, which braces a limb, with electrical stimulation, which restores specific muscle function. Moving beyond bracing, Myo-Orthotics Technology devices restore the functionality of an impaired limb by recreating a natural nerve-to-muscle response using electrical stimulation. WalkAide utilizes advanced sensor technology to analyze leg & foot movement. WalkAide then utilizes Myo-Orthotic Technology to send electrical pulses to the peroneal nerve which prompts the muscles in the leg to raise the foot at the appropriate time.

2007 WalkAide® Training Schedule

DATE	LOCATION
Jul. 13-14	Columbus, OH
Jul. 23-24	Seattle, WA
Aug. 13-14	Dallas, TX
Aug. 27-28	Tampa, FL
Sep. 10-11	Detroit, MI
Oct. 15-16	Newport Beach, C/
Oct. 22-23	Jacksonville, FL
Nov. 5-6	Albany, NY



Experience Our Commitment

Tel: 800.767.7776 x 3 www.spsco.com

Fax: 800.869.7776 estore.spsco.com