



WalkAide is a Class II device, restricted by physician prescription & disbursement through certified practitioners only. SPS is currently registering practitioners for the 2007 WalkAide Training Courses in major cities across the US. Please call your SPS Regional Sales Manager to register for the 2007 WalkAide® Training Courses nearest you, Each course is limited to 25 participants.

WalkAide® System Myo-Orthotics Technology™

Represents the merging of orthotic technology, which braces a limb, with electrical stimulation, which restores specific muscle function. Moving beyond bracing, Myo-Orthotics Technology devices restore the functionality of an impaired limb by recreating a natural nerve-to-muscle response using electrical stimulation. WalkAide utilizes advanced sensor technology to analyze leg & foot movement. WalkAide then utilizes Myo-Orthotic Technology to send electrical pulses to the peroneal nerve which prompts the muscles in the leg to raise the foot at the appropriate time.

2007 WalkAide® Training Schedule

LOCATION	
oort Beach, CA sonville, FL iy, NY	

Come view the New Patient Testimonial Video in Booth #300 at the AOPA National Assembly in Las Vegas.



Tel: 800.767.7776 x 3 www.spsco.com

Fax: 800.869.7776 estore.spsco.com