



Every child  
deserves  
to be a kid!



Model 725

**Compact Double Action Ankle Joint** (Model 725) offers a **variety of options**, all in one complete small package. It comes optimized with pins in the anterior channels and springs in the posterior channels to afford normal gait after initial contact, while providing maximum control and stability of the knee and hip. **Extra springs** and pins are also included to address all the potential needs of a patient's foot/ankle complex.

The Compact Double Action Ankle Joint is supplied **fully assembled** and is designed specifically for use with thermoplastics. It is intended for both small adult and pediatric patients and comes complete with an alignment tool. Designed to provide optimal control in the sagittal plane and increase mediolateral stability, the Compact Double Action Ankle Joint is made with the same **Becker quality** you have come to know and trust since 1933.



p 800-521-2192 | 248-588-7480  
f 800-923-2537 | 248-588-2960

[BeckerOrthopedic.com](http://BeckerOrthopedic.com)   

Available in a modular version:

**Model 726** is a modular version with a larger recess to accommodate 3/16" x 1/2" uprights (not included). The upright attachment holes have also been threaded to allow you to attach your uprights with screws.

