

COYOTE'S TWO CFAB PRODUCTS

BOISE BI-VALVE AFO

- Alternative to the CROW boot and Arizona leather gauntlet brace.
- Increased off-loading of the foot and ankle.
- Fits inside most orthopedic and walking shoes

Good for treatment of a variety of pathologies including Charcot foot, PTTD, osteo arthritis, and protecting ankle fusions post-operatively.



DYNAMIC STRUT AFO

- Some tibial torsion for a more natural gait.
- Responsive energy return.
- Rigid or Medium Flex.

Go to coyotedesign.com/c-fab or call for more information.

OR GET OUR DO IT YOURSELF KIT

AFO Strut Choices

CD205R (Rigid Flex)

Patient Weight Range: 175 - 265lbs

Patient Activity Level: Medium Activity Level

CD205M (Medium Flex)

Patient Weight Range: 125-200lbs

Patient Activity Level: Light to Medium Activity

The strut is designed for average daily use; walking, working, hiking, biking and golf. Has not been tested with running.

Kit is easy to fabricate into a custom AFO



Phone (208) 429-0026
www.coyotedesign.com