Torso-Interface with X-STATIC®

Knit with the magic of silver.



Odor Silver ions X-STATIC work X-STATIC quickly transfers heat from the body and evenly distributes it throughout the garment, conducting the heat to the outside air, and accelerating moisture evaporation.

Solutions that fit.

Building on over 80 years of medical textile innovations, Knit-Rite® brings the most effective antimicrobial agent available and enhanced therapeutic characteristics to the Torso-Interface $^{\circ}$ products with X-STATIC $^{\circ}$ - The Silver Fiber™ from Noble Fiber Technologies.

Like the original seamless body sock, Torso-Interface with X-STATIC may be used with all types of spinal orthoses. Slightly heavier than a T-shirt, it provides better skin protection, conforms to body contours, and also provides the proven benefits of pure silver.

The Warm Weather Solution. Thermally conductive silver fibers in X-STATIC work much like a silver spoon in a cup of hot tea, by conducting heat away from the skin to the ambient air, thereby cooling the body.

X-STATIC also enhances the body's natural ability to evaporate moisture. The conductivity of the X-STATIC fibers acts to accelerate evaporation. As moisture is evaporated, a correspondingly greater amount of moisture can be drawn through the yarn system for subsequent evaporation. Torso-Interface with X-STATIC transports perspiration away from the body faster, creating a more comfortable environment and less potential for the development of dermatitis, bacterial irritation, and odor. As a result, Torso-Interface with X-STATIC allows for increased comfort and better skin health by minimizing moisture contact with the skin and controlling bacteria.

Call Customer Service or your sales representative today for more information on this and our other revolutionary silver products.

Phone: 800-821-3094 Fax: 800-462-4707

Also distributed by:













