

# FREE Rotating Adaptor

## with each purchase of the OFM2 knee

**Leading the world in motion since 2001.**

In celebration of the 10th Anniversary of the OFM2 knee, medi will be giving away our innovative Rotating Adaptor FREE with every OFM2 purchase through June 30, 2011\*. Thank you for continuing to support **medi prosthetics** and providing your patients with the OFM2.



### Our most popular Rotating Adaptor

- Adjusts in half the time: simply loosen any two adjacent alignment screws to rotate, re-tighten and it's fixed.

### OFM2 Knee

- Converts between friction brake and manual locking knee with a simple adjustment.
- Transitions patients between K1-K2.
- Weighs only 495 grams, supports up to 275 pounds.
- Uses a friction brake that releases before unweighting for smoother ambulation.

### medi Cares.

Last year, medi USA sent a team of prosthetists to fit Haitian earthquake victims with OFM2 knees. As a part of ongoing international efforts to rebuild this impoverished nation, medi has donated 500 OFM2 knees to Haitian amputees.

