

*Jacob is more
stable and upright*

JACOB'S PATIENT INSPIRED SOLUTION

Ultraflex addresses the multiple treatment goals involved with crouch. Jacob wears an Ultraflex dynamic assist stretching KO (at rest/nighttime) to maintain and increase ROM. Ultraflex Adjustable Dynamic Response™ (ADR™) AFOs (daytime) are worn bilaterally to control dorsiflexion, provide normal range, increase stability, and improve balance. The Ultraflex stretching KO worn at night/at rest is easily connected/disconnected to and from the ADR™ AFO with the Ultraflex UltraQuick Release™. The stretching KO/AFO ADR™ combination is billable as one brace. The family reports that Jacob is getting the stretching he needs, has improved balance and stability, and is more upright—all of which give Jacob greater function.



For education for you and your referrals, please call:
800-220-6670

Ultraflex[®]
Patient Inspired Solutions
www.ultraflexsystems.com

Existing coding and coverage applies to all ADR™ technology.

Thank you to Keith Smith, CO, LO, FAAOP of O&P Labs, St. Louis, MO and Courtney Dunn, PT, DPT, of St. Louis Children's Hospital for their clinical contributions. A special thanks to Jacob and his family.